# **Specimen Collection and Submission Guidance for**

Chlamydia trachomatis and Neisseria gonorrhoeae Screening (by Amplified RNA Probe)

Collect each specimen using the appropriate collection and transport kit. Follow the manufacturer's (Hologic®) specimen collection instructions.

#### Vaginal, Rectal, and Throat Swabs

Use: Aptima® Multitest Swab Specimen Collection and Transport Kit

#### **Specimen Storage and Shipping**

**Vaginal Swabs:** After collection, store and ship swabs in specimen transport tube at 2°C to 30°C.

**Rectal and Throat Swabs:** After collection, store and ship swabs in specimen transport tube at 4°C to 30°C.

- Swabs must be received within 58 days of collection.
- Freeze vaginal swabs at -20 °C to -70 °C within 7 days of collection for long-term storage.

# Actions 144

Aptima® Multitest Swab Specimen Collection and Transport Kit

# **Urine Specimens**

**Use:** Aptima® Urine Specimen Collection and Transport Kit

#### **Specimen Storage and Shipping**

Transfer urine sample to Aptima<sup>®</sup> urine specimen transport tube within 24 hours of collection.

Store and ship transport tubes at 2°C to 30°C.

- Specimens must be received within 28 days of collection.
- Freeze at -20 °C to -70 °C within 7 days of collection for longterm storage.



Aptima® Urine Collection and Transport Kit

### **Endocervical and Urethral Swabs**

Use: Aptima® Unisex Swab Specimen Collection and Transport Kit

#### **Specimen Storage and Shipping**

After collection, store and ship swabs in specimen transport tube at 2°C to 30°C.

- Swabs must be received within 58 days of collection.
- Freeze at -20 °C to -70 °C within 7 days of collection for longterm storage.



Aptima® Unisex Swab Specimen Collection and Transport Kit

# **Specimen Collection and Submission Guidance for**

Testing of C. trachomatis (CT) and N. gonorrhoeae (GC) by Amplified RNA Probe

#### **Ensure Specimen Labels Have at Least Two Unique Identifiers**

Three unique patient identifiers on specimen are preferred.




Three patient identifiers provided on this label.

- Name
- 2. Date of Birth
- 3. **Medical Record Number**

#### Provide Patient Identifiers in Sections 2 and 3 of Form G-2B

Patient identifiers on specimen label and G-2B form must match exactly.

**Date of Collection must** be provided in Section 3.



# **Identify Specimen Type in Section 3**

Select only **one** specimen source.

#### Questions About . . .

Specimen Collection/Suitability:

Specimen Shipping:

**Supply Ordering:** 

(512) 776-7657 or 512-776-2449

(512) 776-7598 or 1-888-963-7111 ext. 7578 (toll free)

(512) 776-7661 or ContainerPrepGroup@dshs.texas.gov

# Submitter Accounts, Submission Forms, or Result Reports: (512) 776-7578 or LabInfo@dshs.texas.gov

# **Request Test in Section 4.2**

Select "GC/CT, amplified RNA probe".

	4.2 Bacteriology	
_	Clinical Specimen  ☐ Aerobic Isolation ☐ Anaerobic Isolation	Definitive Idea  ☐ Anaerobic ide  Organism Su
	☐ Culture, stool ☐ Diphtheria Screen ▼ GC/CT, amplified RNA	☐ Bacillus spp.
	probe	☐ Campylobacti ☐ Enteric bacter ☐ Gram Negativ

### **Select the Payor in Section 6**

Check the **appropriate** box as Payor. If left empty, the submitter is charged.

RED	☐ Medicaid (2)  Medicaid/Medicare #:	☐ Medicare (8)
** REQUIRED	☐ Submitter (3) ☐ BIDS (1720) ☐ BT Grant (17[9) ☐ HIV / STD (1608) ☐ IDEAS (1610)	☐ Immunizations (1609) ☐ Private Insurance* (4) ☐ TIPP (5144) ☐ Zoonosis (1620) ☐ Other: