Specimen Collection and Submission Guidance for

Norovirus PCR Test

Molecular Tracking of Norovirus to Better Understand Outbreaks and Transmission

Norovirus Stool Specimen Collection

Required Specimen:

- Raw, fresh stool
- Raw, frozen stool
- Stool in Cary-Blair transport medium

Required Volume:

- **Raw Stool:** Minimum 400 μL, but more is preferred.
 - Use a sterile container with tight-fitting or screw-top lid.
- Cary-Blair Stool: Follow manufacturer's guidelines.

Required Storage and Shipping Temperature:

- Raw Stool: Ship frozen at -20°C on dry ice, or cold at 2°C–8°C on ice packs.
- Cary-Blair: May be shipped cold or at ambient temperatures.

Norovirus Shipping and Labeling Requirements

Ship as: Category B Biological Substance, UN3373 Specimen must be:

- Triple packaged to withstand shock, pressure changes, leaks, and other ordinary handling conditions while in transit.
- Packaged with enough absorbent material such as cellulose wadding to soak up the contents of the specimen container.
- **Shipped overnight** in sealed, insulated containers with cold packs (refrigerated specimens), or dry ice (frozen specimens).
- Ensure all containers are securely closed to prevent leaks.
 - Secure lids by wrapping in paraffin film (e.g., Parafilm).
- Ensure cold specimens are packed with multiple ice packs.
- Ensure transport media is not expired.

Visit DSHS' online <u>Specimen Shipping and Mailing Guidance</u> for more details on shipping Category B substances to the Laboratory. <u>Specimens received out of temperature range will be rejected.</u>



For specimen size comparison, a single garden pea is approximately 200 μL in volume. Image source: pixabay.com





Specimen Collection and Submission Guidance for

Norovirus PCR Test

Norovirus Specimens Must be Labeled and Submitted with a G-2B Submission Form

Label Specimen With Unique Identifiers

Every specimen must have at least two unique patient identifiers on its label.




3 06161858

- Three patient identifiers provided on this label:
- 1. Name
- 2. Date of Birth
- 3. Medical Record Number

Provide Patient Identifiers in Sections 2 and 3 of Form G-2B

Patient identifiers on specimen label and G-2B submission form must match.

Date of Collection must be provided in Section 3.



Select Specimen Type in Section 3

Check "Feces/ stool" as the

specimen type.

☐ Eve Swah ☐ Feces / stool ☐ Gastric (Aspirate)

Select Test Type in Section 4.4

Check "Norovirus" under *Molecular Studies*.

Select "IDEAS (1610)" in Section 6

Check the IDEAS (1610) box as the Payor.

-			
I		☐ Medicaid (2)	☐ Medicare (8)
ı	۵	Medicaid/Medicare #:	
ı	REQUIRED	_	
ı	5	☐ Submitter (3)	☐ Immunizations (1609)
ı	g	☐ BIDS (1720)	☐ Private Insurance* (4)
ı		☐ BT Grant (1719)	☐ TIPP (5144)
ı	*	☐ HIV / STD (1608)	☐ Zoonosis (1620)
ı		IDEAS (1610)	☐ Other:
L			

NOTE: If the 'Date of Collection' field is not completed, the specimen will be rejected. Date of Collection (mm/dd/yyyy) ** Time of Collection *** Collection ** Collection ** PM Unique Identification Number ** Comments or Additional ID e.g., MRN / Alien # / Accession ID 06161858

Questions About...

Specimen Collection/Suitability: (512) 776-6510

Specimen Shipping: (512) 776-7598 or 1-888-963-7111 ext. 7578 (toll free)

Norovirus Surveillance Program: <u>FoodborneTexas@dshs.texas.gov</u>

Submitter Accounts, Submission Forms, or Result Reports: (512) 776-7578 or LabInfo@dshs.texas.gov