# Specimen Collection and Submission Guidance for Measles (Rubeola) Serum Specimens

#### Submitting Serum Specimens for Detection of IgG/IgM Antibodies to Measles

#### **Measles Serum Specimen Collection and Storage**

**Report suspected measles cases immediately.** Reporting contacts by county/public health region available <u>here</u>.

Required Specimen: Serum

Required Volume: 200 µL, but more is preferred

#### **Required Storage and Shipping Temperatures:**

- Store and ship cold at 2°C-8°C if specimen will arrive at the Laboratory within 48 hours of collection.
- Store and ship frozen at -20°C or colder if specimen will arrive at Laboratory more than 48 hours after collection.

Ensure specimen collection kits are not expired!

#### **Specimen Shipping and Labeling Requirements**

**Ship as:** Category B Biological Substance, UN3373 Specimens must be:

- **Triple packaged** to withstand shock, pressure changes, leaks, and other ordinary handling conditions while in transit.
- **Packaged with enough absorbent material** in secondary container to soak up the contents of the primary container.
- **Shipped overnight** in insulated containers with frozen cold packs (for cold specimens), or dry ice (for frozen specimens).

Ensure containers are securely closed to prevent leaks.

- Secure specimen lids by wrapping in paraffin film (e.g., Parafilm).
- **Pack** cold specimens with multiple ice packs.

Visit DSHS' online <u>Specimen Shipping and Mailing Guidance</u> for more details on shipping Category B substances.

#### Ensure "Saturday Delivery" box is selected for Friday shipments. Do not ship for Sunday or holiday delivery.

Specimens received out of temperature range will be rejected.



A child's cheek showing the characteristic measles rash. (CDC/ Tatiana Lanzieri, MD, MPH, 2024)



**NEVER FREEZE Serum** 

in Separator

**Tubes!** 

Ship Measles Specimens Cold or Frozen Overnight!



lealth and Human Texa ervices Hea

TXAS

Texas Department of State Health Services

dshs.texas.gov

Updated 02/2025

# Specimen Collection and Submission Guidance for Measles (Rubeola) Swab Specimens

Submitting Measles Swab Specimens for Viral RNA Detection

#### Measles Swab Specimen Collection and Storage

# **Report suspected measles cases immediately.** Reporting contacts by county/public health region available here.

#### **Required Specimens:**

- Preferred: Throat swabs
- Acceptable: Nasopharyngeal or oral swabs
   Use dacron or rayon swabs only; no cotton tips or wooden shafts.
- Use Universal Transport Media or Viral Transport Media.

#### Required Volume: 1 mL-2 mL

#### **Required Storage and Shipping Temperature:**

- Store and ship cold at 2°C—8°C if specimen will arrive at Laboratory within 48 hours of collection.
- Store and ship frozen at -70°C if specimen will arrive at Laboratory more than 48 hours after collection.

Ensure specimen collection kits are not expired!

#### **Specimen Shipping and Labeling Requirements**

**Ship as:** Category B Biological Substance, UN3373 Specimen must be:

- Triple packaged to withstand shock, pressure changes, leaks, and other ordinary handling conditions while in transit.
- **Packaged with enough absorbent material** in secondary container to soak up the contents of the primary container.
- **Shipped overnight** in insulated containers with enough cold packs to keep them cold or enough dry ice to keep them frozen.
  - Pack cold specimens with multiple frozen ice packs.
- Ensure containers are securely closed to prevent leaks.
  - Secure lids by wrapping in paraffin film (e.g., Parafilm).


Visit DSHS' online <u>Specimen Shipping and Mailing Guidance</u> for more details on shipping Category B substances.

Ensure "Saturday Delivery" box is selected for Friday shipments. Do not ship for Sunday or holiday delivery. Specimens received out of temperature range will be rejected.



EXAS

Texas Department of State Health Services



(CDC/ Allison M. Maiuri, MPH, CHES (2016))

Ship Measles Specimens Cold or Frozen Overnight!

dshs.texas.gov

Updated 02/2025


# **Specimen Collection and Submission Guidance for**

**Measles Serum Specimens for Antibody Detection** 

\*\*\*All Facilities Must Have a DSHS Submitter Account to Submit Specimens\*\*\*

### Label Specimen With Unique Identifiers

Every specimen must have at least two unique patient identifiers on its label.




3. Medical Record Number

## Provide Patient Identifiers in Sections 2 and 3 of Form G-2A

**Health Services** 

Services

Patient identifiers on specimen label and G-2A submission form must match. **Date of Collection must** be provided in Section 3.



Updated 02/2025

# **Specimen Collection and Submission Guidance for**

**Measles Swab Specimens for Viral RNA Detection** 

\*\*\*All Facilities Must Have a DSHS Submitter Account to Submit Specimens\*\*\*


## Label Specimen With Unique Identifiers

Every specimen must have at least two unique patient identifiers on its label.



#### Provide Patient Identifiers in Sections 2 and 3 of Form G-2V

Patient identifiers on specimen label and G-2V submission form must match. Date of Collection must be provided in Section 3.



#### Questions About ...

#### Laboratory Mailing Address:

Public Health Laboratory Division, MC 1947 Department of State Health Services 1100 W. 49th Street, Austin, TX 78756-3199

Specimen Collection/Suitability: Measles Surveillance Program: Specimen Shipping: Submitter Accounts, Submission 512-776-7594 or 512-776-7760 <u>Measles (Rubeola) | Texas DSHS</u> 512-776-7598 or 1-888-963-7111 ext. 7578 (toll free)

Submitter Accounts, Submission Forms, or Result Reports: 512-776-7578 or LabInfo@dshs.texas.gov



Texas Department of State Health Services


## **Select Specimen Type in Section 3**

<b>Check</b> the appropriate specimen type.	Feces/stool     Nasopharyngeal swab     Nasal Swab	
Select Test Type in Section 4		
Check		

CHECK	
Measles PCR.	Measles PCR
Identify	Vaccine Received:  Yes  No Date Vaccine Received:
vaccine status,	Travel History (if known):
travel history.	

# **Select Payor in Section 6**

**Check** the appropriate Payor. Do not leave empty.



dshs.texas.gov