



TEXAS
Health and Human
Services

Texas Department of State
Health Services

Texas Department of State Health Services
Tuberculosis and Hansen's Disease Unit

Tuberculosis Specimen Shipping Guide

Created September 1, 2020; Revised October 24, 2024

Table of Contents

FedEx Accounts	3
Shipping TB Specimen to Designated Laboratories	4
DSHS Laboratory	4
DSHS South Texas Laboratory (STL)	4
University of Florida – Infectious Disease Pharmacokinetics Laboratory (IDPL)...	5
Quest Diagnostics	5
Shipping Biological B Specimen and Supplies.....	6
Common Reasons for Unsatisfactory Specimen	8
Shipping Labels	9
Category B – Biological Substances Label UN 3373.....	9
Dry Ice – Class 9 Miscellaneous Dangerous Goods UN 1845.....	9
Frequently Asked Questions for Shipping to DSHS Laboratories	10
G-MYCO Specimen Submission Form – DSHS Austin	12
G-2A Serology Specimen Submission Form- DSHS Austin	13
F40-TB Elimination Specimen Submission Form - DSHS South Texas Laboratory (STL)	14
F40-B Specimen Submission Form - DSHS South Texas Laboratory (STL).....	15

FedEx Accounts

Account Set-Up

Texas Department of State Health Services (DSHS) public health region (PHR) and DSHS-contracted local health department (LHD) tuberculosis (TB) programs must establish a courier account with DSHS TB and Hansen's Disease Unit (TB Unit) to ship TB Unit-approved specimen to DSHS laboratories and other DSHS-contracted laboratories using Federal Express (FedEx).

To set up a new account, regional and local health department (R/LHD) TB programs may notify the TB Unit's Administrative Team at TBProgram@dshs.texas.gov by submitting the [Notice of Change in TB Personnel](#) form with the following information:

- Name of submitter
- Work Address (including zip code)
- Email address of clinic contact
- Name of public health region/county/clinic

An account must be established for each DSHS-contracted LHD clinic site and PHR field office. Once established, programs may use this account to ship all TB specimen to the laboratories outlined in this guide.

Courier Services Offered

1. FedEx Priority Overnight: *Select for specimen shipping.*
2. FedEx Standard Overnight: *Select for return labels.*
3. FedEx Home Delivery
4. FedEx Ground

Helpful FedEx Resources

- Customer Service Telephone Number: 1(866) 477-7529
- How to Ship Clinical Samples (prevent leaking specimen):
 - <https://www.fedex.com/en-us/shipping/how-to-ship-clinical-samples.html>
- Service Guide:
 - https://www.fedex.com/content/dam/fedex/us-united-states/services/Service_Guide_2024.pdf
- Ground Service Map (transit):
 - <http://www.fedex.com/grd/maps/ShowMapEntry.do>
- InSight (advanced shipment tracking):
 - <https://www.fedex.com/en-us/tracking/insight.html>
- Return Shipments (labels):
 - <https://www.fedex.com/en-us/service-guide/return-shipments.html>

Shipping TB Specimen to Designated Laboratories

R/LHD TB programs may order free specimen shipping boxes from the TB Unit. These boxes may not be used to ship specimen for other programs. TB specimen may be shipped to the laboratories outlined in this document.

DSHS Laboratory

1100 West 49th Street Austin, Texas 78756

Phone: (512) 776-7318 or (512) 776-7598; Fax (512) 776-7294

<https://www.dshs.texas.gov/lab/default.shtm>

- **Tests performed:**
 - Acid Fast Bacilli (AFB) smear and culture
 - Nucleic Acid Amplification Test (NAAT)
 - Drug susceptibility tests (DSTs)
 - HIV, Hepatitis B and C
- Ship all specimen by FedEx Priority **Overnight** Monday through Thursday (so shipment does not arrive on a weekend).
- Do not ship on Friday, Saturday, or the day prior to a holiday.
- **DSHS laboratory specimen collection criteria:**
<https://www.dshs.texas.gov/laboratory-services/programs-laboratories/microbiology-unit/mycobacteriology-mycology/growth-detection-mycobacteria>
- **DSHS Standing Delegation Orders (SDOs) specimen collection procedures:**
 - [TB Sputum Collection SDO](#)
 - [TB Blood Specimen SDO](#)
- **Ordering forms for sputum collection supplies:**
<https://www.dshs.texas.gov/laboratory-services/laboratory-testing-services-manual-forms-laboratory-fee-schedule>

DSHS South Texas Laboratory (STL)

1301 S. Rangerville Road, Harlingen, TX 78552

Phone: (956) 364-8746 or (956) 364-8753 (TB); (956) 364-8751 (Hematology); and (956) 364-8752 (Clinical); Fax (956) 412-8794

https://www.dshs.texas.gov/lab/so_tx_lab.shtm

- **Tests performed:**
 - Blood testing results for chemistry, special chemistry, hematology
 - *Mycobacteriology testing (operations are on hold until further notice):*
 - *AFB smear and culture (include urine testing), NAAT, and DSTs*
- Ship all specimen by FedEx Priority **Overnight** Monday through Thursday.
- Do not ship on Friday, Saturday or the day prior to a holiday.

University of Florida – Infectious Disease Pharmacokinetics Laboratory (IDPL)

Infectious Disease Pharmacokinetics Laboratory

Dr. Charles Peloquin, Pharm D.

University of Florida

1600 SW Archer RD., P4-30

Gainesville, FL 32610

Phone: (352) 273-6710

<https://idpl.pharmacy.ufl.edu/>

- **Test performed:**
 - Therapeutic drug monitoring
- Pack samples upright in Styrofoam box and ship FedEx Priority **Overnight**.
- Only ship specimen Monday through Wednesday. Do *not* ship on Thursday through Sunday, or the day prior to a holiday.
- Package properly for dry ice handling, including a dry-ice specific label. Ship on at least five pounds of dry ice.
- Include a FedEx return label inside the cold box in a plastic bag to prevent damage or loss during transit.
- **Shipping details:**
<https://idpl.pharmacy.ufl.edu/wordpress/files/2022/08/Instructions-sample-handling-UFSHands-v-08.22.pdf>.
- **DSHS TB Unit's Therapeutic Drug Monitoring Process:**
<https://www.dshs.texas.gov/sites/default/files/IDCU/disease/tb/forms/PDFS/TherapeuticDrugMonitoringProcess.pdf>.

Quest Diagnostics

Quest Client Services: 866-MYQUEST (866-697-8378)

www.questdiagnostics.com


- **Tests performed:**
 - Interferon Gamma Release Assays (IGRAs):
 - T-SPOT®-TB Test
 - QuantiFERON®-TB Gold Plus One Tube Test (QFT)
- R/LHDs can request a Quest courier pick up either via Quantum (the "Specimen Pick Up" button) or by calling 1-866-MY QUEST.
- Both IGRAs must be picked up by the Quest courier or shipped with Quest's FedEx label and packaging the same day they are drawn.
- **Quest T-SPOT details:**
 - <https://testdirectory.questdiagnostics.com/test/test-detail/37737/t-spottb?cc=MASTER>
- **Quest QFT details:**
 - <https://www.questdiagnostics.com/healthcare-professionals/clinical-education-center/faq/faq204>

Shipping Biological B Specimen and Supplies

Sputum and blood are categorized as “Biological Substance, Category B” for shipping. This means they are infectious substances transported for diagnostic purposes and submitters must adhere to shipping requirements. Details are at https://www.fedex.com/content/dam/fedex/us-united-states/services/UN3373_fxcom.pdf.

TB Programs must ship Category B specimen with three packing layers (see Figure 1). Collect specimen in a **primary receptacle** (such as vacutainer for blood or sterile blue top tube for sputum), wrap in absorbent packaging (such as tissue or cotton), place in a **secondary receptacle** (a leak proof container), and ship in a rigid outer covering (a box, or **third receptacle**) with frozen gel packs for cold shipping or with dry ice for therapeutic drug monitoring. Specimen require labels for Category B and dry ice, when applicable.

Figure 1: Packing Category B Specimen for Transport



Source: <https://www.cdc.gov/smallpox/lab-personnel/specimen-collection/pack-transport.html>.
See https://www.dshs.texas.gov/lab/mrs_shipping.shtm#Samples for DSHS shipping recommendations.

Primary Receptacle for Sputum



Primary Receptacle for Blood



Secondary Receptacle(s): Hard Plastic Liner Tube or 95 KPa Biohazard Bag



Third Receptacle: Rigid Outer Covering



Gel Packs, Kept Frozen for Cold Shipping



Biological B Specimen Shipping Label



Dry Ice Shipping Label



Follow these shipping steps for each specimen:

1. Wrap the primary receptacle in absorbent material (i.e. a paper towel or cotton) and place into the secondary receptacle. Ensure there is enough absorbent material around and on top of the primary receptacle tube so that it cannot easily move around.
2. Include the completed laboratory requisition with the specimen and ensure that it will not get wet. Do not put patient information on outer or secondary container or lids.

3. Place the secondary receptacle with the enclosed/affixed requisition in the third receptacle that is made of a Styrofoam layer and an outer cardboard box (rigid outer covering). Place at least two icepacks in the box, including one on bottom and one on top to “sandwich” the specimen, or at least five pounds of dry ice if sending therapeutic drug monitoring specimen. Place absorbent paper towels on the ice to ensure melting does not wet the bag or laboratory requisition.
4. Place the FedEx label on the inside flap in a pouch to prevent damage or loss during transit or when the box is opened. If the return label is placed on top of the Styrofoam container (inside the taped seam), it may be sliced in half when staff open the box.
5. Close the box and tape securely unless otherwise directed by the courier. Affix the “To” address label and **UN 3373 Biological B Specimen label** (and **UN 1845 Dry Ice Label, when applicable**) to the outside of the box where clearly visible.

Common Reasons for Unsatisfactory Specimen

- 1. Leaking specimen.** This occurs when the blue lid for sputum samples is not fully shut, causing the specimen to leak during transit. To avoid leaking specimen, tighten the blue lid carefully. Ensure the lid is threaded correctly and tightened; do not overtighten or strip the threads. Pack the blue-top tube tightly inside the black top outer receptacle with absorbent material to reduce the risk of leakage.
- 2. Missing patient identifiers.** Labs will reject specimen without two matching patient identifiers on the requisition form and on the outside of the primary receptacle. Common identifiers are **patient name (first and last)** and **date of birth**, or **patient name** and **medical record number**.


For other examples of unsatisfactory specimen sent to South Texas Laboratories (STL), see <https://www.dshs.texas.gov/laboratory-services/south-texas-laboratory/stl-unsatisfactory-specimen-examples>.

Shipping Labels

Category B – Biological Substances Label UN 3373

“Biological Substance, Category B” must appear in 6-mm-high text on the outer package adjacent to a diamond-shaped mark. The UN 3373 mark must be in the form of a square set at an angle of 45 degrees. Each side of the UN 3373 diamond should measure a minimum of 2" (50 mm). The width of the diamond rule line must be a minimum of 2 mm, and the letters and numbers must be at least 6 mm high.”

These labels may be printed or purchased from various online sources. Contact FedEx for labeling support at https://www.fedex.com/content/dam/fedex/us-united-states/services/UN3373_fxcom.pdf.



Dry Ice – Class 9 Miscellaneous Dangerous Goods UN 1845

An International Air Transport Association (IATA) Class 9 Miscellaneous label must appear on all dry ice shipments. FedEx offers a dry ice label that when correctly completed, satisfies the IATA marking and labeling requirements.

The following permanent markings are required on the outer packaging of all IATA dry ice shipments:

- Dry Ice
- UN 1845
- Net weight of dry ice in kilograms
- Name and address of the shipper
- Name and address of the recipient



Print this label at https://www.fedex.com/content/dam/fedex/us-united-states/services/Dry_Ice_Label.pdf.

Frequently Asked Questions for Shipping to DSHS Laboratories

What FedEx shipping supplies are needed, and where can R/LHDs find them?

After setting up the R/LHD FedEx account, shipping labels can be printed locally. R/LHDs can also order clear plastic envelopes (to adhere shipping label to the cardboard cold box) and preprinted shipping labels from FedEx:

<https://www.fedex.com/en-us/shipping/packing.html>.

Where do R/LHDs order cold boxes and gel ice packs?

Insulated cardboard cold boxes and gel ice packs that can be frozen are available via the DSHS Pharmacy Inventory Ordering System (PIOS). When searching for bulk options, look for "Specimen Transport Cold Box" and "Specimen Transport Gel Packs".

Should cold boxes be reused?

Yes, please reuse cold boxes. Include a FedEx return label with *all* shipments to ensure the box is returned to the sending R/LHD program. The return label should include the account information provided by the TB Unit. Ensure the return address is accurate. If submitters are unsure of their account details, please email TBProgram@dshs.texas.gov. Maintain boxes in good condition by preventing specimen leaks and tears.

Where can R/LHDs order specimen collection supplies?

Sterile TB specimen collection tubes, secondary containers, and cardboard outer mailers are available via the DSHS laboratory. Download the most current order form from: <https://www.dshs.texas.gov/laboratory-services/laboratory-testing-services-manual-forms-laboratory-fee-schedule>.

What type of TB specimen should be sent in the cold boxes?

Cold boxes are used to ship TB specimen to DSHS laboratories in Austin or South Texas, or to outside reference laboratories for TB testing. This includes sputum samples and blood tests. The boxes are insulated for shipment of cold specimen when used with gel ice packs, or room temperature specimens when gel ice packs are not needed. See https://www.dshs.texas.gov/lab/MRS_specimens.shtm for details.

Note: Do not use the cold box program for shipping specimen to Quest. They provide shipping supplies including pick-up service to ship blood specimen for IGRA testing to their laboratories.

How many specimens fit in one cold box?

Sites may ship 50 ml of specimen or fewer per box. Two or more primary receptacles may be included per box.

How many cold packs are recommended per sputum canister or box?

Typically, two gel ice packs per box is recommended. However, FedEx does not use temperature-controlled trucks to transport boxes, so consider outside temperatures and transport drive time. For extreme heat, use three or four gel ice packs.

Do R/LHDs need to ship all sputum to the DSHS laboratory in a cold box via FedEx, or can R/LHDs send samples via regular mail in the brown mailing canisters?

Now that the TB Unit provides a cold-box FedEx account for all TB programs, it is recommended that every sputum sample is collected and shipped via cold box with gel ice packs. This will provide DSHS laboratories the best possible specimen to test, as it will arrive quickly via FedEx at the recommended cold temperature necessary for testing.


There may be times when this is not possible. For example, cold-box shipping is not possible when mailing canisters are left for the patient to self-collect and send via U.S. mail. When this occurs, instruct patients to keep the sample refrigerated before shipping in the brown outer mailer. Self-collected specimens should be saved, refrigerated and shipped in the cold boxes with ice packs when public health personnel can pick up the specimen from the patient.

How often does DSHS return shipping boxes? How can R/LHDs ensure that DSHS returns boxes?

DSHS returns cold boxes on Tuesday through Friday afternoons, 1-2 days after they have been picked up from FedEx. For example, if the submitter sends the cold box on Monday, FedEx will deliver the box on Tuesday morning and DSHS will return it to the submitter on Tuesday afternoon (or Wednesday, depending on distance). The only exceptions are closed holidays which would extend the return to the next business day.

Cold boxes received by DSHS laboratories are returned in brown cardboard boxes to prevent wear and tear on the outer cardboard box that protects the Styrofoam containers, provided there is a return shipping label. If no return label is received with the cold box and no information is available inside the cold box to determine who it belongs to, the lab will not return the cold box. Writing "please return" is not acceptable. The return label should include the name and address of facility, phone number and contact name.

G-MYCO Specimen Submission Form – DSHS Austin

 TEXAS Department of State Health Services		G-MYCO Specimen Submission Form (Jan 2020) CAP# 3024401 CLIA 45D0660644 www.dshs.texas.gov/lab		*** FOR DSHS USE ONLY ***	
Specimen Acquisition: (512) 776-7598 Section 1. SUBMITTER INFORMATION (** REQUIRED, DO NOT ALTER) Submitter/FPI Number ** Submitter Name ** NPI Number ** Address ** City ** State ** Zip Code ** Phone ** Contact Fax ** Clinic Code				Section 6. ORDERING PHYSICIAN INFORMATION – (** REQUIRED) Ordering Physician's Name ** Ordering Physician's License # **	
Section 2. PATIENT INFORMATION – (** REQUIRED) NOTE: Patient name MUST match name on this form, Medicare/Medicaid card & specimen container. Specimen must have two (2) identifiers that match this form. Last Name ** First Name ** MI Address ** Residence Address City ** State ** Zip Code ** Country of Origin / Bi-National ID # DOB (mm/dd/yyyy) ** Sex Race: <input type="checkbox"/> White <input type="checkbox"/> Black or African American <input type="checkbox"/> American Indian / Native Alaskan <input type="checkbox"/> Asian <input type="checkbox"/> Native Hawaiian / Pacific Islander <input type="checkbox"/> Other Ethnicity: <input type="checkbox"/> Hispanic <input type="checkbox"/> Non-Hispanic <input type="checkbox"/> Unknown Pregnant? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Unknown Date of Collection ** (REQUIRED) Time of Collection: <input type="checkbox"/> AM <input type="checkbox"/> PM Medical Record # ICD Diagnosis Code (1) ICD Diagnosis Code (2) ICD Diagnosis Code (3) Collected by:				Section 7. FAVOR SOURCE – (REQUIRED) 1. Reflex testing will be performed when necessary and the appropriate party will be billed. 2. If the patient does not meet program eligibility requirements for the test requested and no third party payor will cover the testing, the submitter will be billed. 3. Medicare generally does not pay for screening tests. Please refer to applicable Third Party payor guidelines for instructions regarding covered tests, benefit limitations, medical necessity determinations and required beneficiary Notice (ASN) requirements. 4. If Medicaid or Medicare is indicated, the appropriate member ID number is required. Please write it in the space provided. 5. If private insurance is indicated, the responsible party information below is designated with an asterisk (*). 6. Check only one box below to indicate who should bill the submitter, Medicaid, Medicare, private insurance, or TB Program. <input type="checkbox"/> Medicaid (2) <input type="checkbox"/> Medicare (5) Medicaid/Medicare #: <input type="checkbox"/> Submitter (3) <input type="checkbox"/> Private Insurance (4) <input type="checkbox"/> Other (1) (1610) <input type="checkbox"/> IDEAS (1610) <input type="checkbox"/> TB Elimination (1619) <input type="checkbox"/> Other: MO/Insurance Care / Insurance Company Name * Responsible Party * Insurance Phone Number * Responsible Party's Insurance ID Number * Group Name * Group Number * *I hereby authorize the release of information related to the services described here and hereby assign any benefits to which I am entitled to the Texas Department of State Health Services, Laboratory Services Section. Signature of patient or responsible party. Signature * Date *	
Section 3. SPECIMEN SOURCE OR TYPE – (** REQUIRED) <input type="checkbox"/> Abdominal fluid <input type="checkbox"/> Eye <input type="checkbox"/> Sputum: Natural <input type="checkbox"/> Ascites (site) <input type="checkbox"/> Feces/Stool <input type="checkbox"/> Fluorescent <input type="checkbox"/> Aspire (site) <input type="checkbox"/> Gastric <input type="checkbox"/> Tissue (site) <input type="checkbox"/> BAL <input type="checkbox"/> Lesion (site) <input type="checkbox"/> Urinal <input type="checkbox"/> Biopsy (site) <input type="checkbox"/> Lymph node (site) <input type="checkbox"/> Urine (site) <input type="checkbox"/> Bronchial washings <input type="checkbox"/> Nasopharyngeal <input type="checkbox"/> Other <input type="checkbox"/> Cervical <input type="checkbox"/> Pleural Fluid/PLF <input type="checkbox"/> CSF <input type="checkbox"/> Sputum: Induced				Section 8. SUSCEPTIBILITY TESTING Is MDR M. tuberculosis suspected? <input type="checkbox"/> Yes <input type="checkbox"/> No Note: Drug susceptibility tests are performed automatically on patient's initial M. tuberculosis isolate. <input type="checkbox"/> MTB Primary Drug Susceptibility Panel Plus Ofloxacin: <input type="checkbox"/> Ethambutol <input type="checkbox"/> Isoniazid <input type="checkbox"/> Pyrazinamide (PZA) <input type="checkbox"/> Rifampin <input type="checkbox"/> Ofloxacin <input type="checkbox"/> MTB PZA Susceptibility Test Only: <input type="checkbox"/> Capreomycin <input type="checkbox"/> Ethionamide <input type="checkbox"/> Isoniazid <input type="checkbox"/> Kanamycin <input type="checkbox"/> Clofazimine <input type="checkbox"/> Rifabutin <input type="checkbox"/> Rifampin <input type="checkbox"/> Clofazimine <input type="checkbox"/> M. kansasii Susceptibility Test: <input type="checkbox"/> Agar, Rifampin	
Section 4. CLINICAL SIGNIFICANCE FOR RAW UNPROCESSED SPECIMENS: <input type="checkbox"/> AFB Smear Only (for release from isolation) <input type="checkbox"/> AFB Smear and Culture <input type="checkbox"/> AFB Smear, Culture and Direct Nucleic Acid Amplification (Laboratory Diagnostic Specimens Only) FOR PROCESSED SPECIMENS ONLY: <input type="checkbox"/> For Respiratory Diagnostic Specimens <input type="checkbox"/> Direct NAAT (M. tuberculosis) NAAT ONLY – NO CULTURE PERFORMED Please provide the AFB Line result for this processed specimen: <input type="checkbox"/> For AFB Smear Positive Specimen <input type="checkbox"/> Direct HPLC for Mycobacterium species, not M. tuberculosis **** Prior authorization required **** Telephone (512) 776-7542 for authorization.				Section 5. REFERRED PURE CULTURE <input type="checkbox"/> Referred AFB isolate identification <input type="checkbox"/> MTB Genotyping Only/for Compliance <input type="checkbox"/> Fungal isolate identification <input type="checkbox"/> Actinomyces, Aerobic, identification NOTE: Please see the form's instructions for details on how to complete this form. Visit our web site at www.dshs.texas.gov/lab . All dates must be entered in mm/dd/yyyy format.	
FOR LABORATORY USE ONLY				Specimen Received: <input type="checkbox"/> Room Temp. <input type="checkbox"/> Cold <input type="checkbox"/> Frozen Laboratory Services Section: 1100 West 49th St Austin, TX 78756	

- Use the **G-MYCO Specimen Submission Form** for mycobacteriology and TB specimen testing.
- **Section 1.** Ensure all information is updated and current.
- **Section 2.** Fill out completely. Patient name must match exactly the patient name labeled on the specimen container.
- **Section 3.** Specimen source must be provided. This will direct how the specimen is processed in the laboratory.
- **Section 4.** Requested test must be selected or specimen is unsatisfactory for testing.
- **Sections 5 and 8.** These sections are only for referred isolates from other laboratories.
- **Sections 6 and 7.** Fill out completely. Ensure *TB Elimination* is marked in "Payor Source" in Section 7.

Visit laboratory website for the most recent requisition version
https://dshs.texas.gov/lab/MRS_forms.shtm

G-2A Serology Specimen Submission Form- DSHS Austin

G-2A Specimen Submission Form (June 2020) Rev. 1		****For DSHS Use Only****	
TEXAS Department of State Health Services Specimen Acquisition: (512) 778-7598		CAP# 302440 CLIA #45D0690644 www.dshs.texas.gov/lab	
Section 1. SUBMITTER INFORMATION (** REQUIRED)		Section 7. ORDERING PHYSICIAN INFORMATION (** REQUIRED)	
Submitter/TPN Number **		Ordering Physician's NPI Number **	
Submitter Name **		Ordering Physician's Name **	
NPI Number **			
Address **			
City **	State **	Zip Code **	
Phone **	Contact		
Fax **	Clinic Code		
Section 2. PATIENT INFORMATION (** REQUIRED)			
NOTE: Patient name on specimen MUST match name on this form & Medicaid/Medicare card. Specimen must have two (2) identifiers that match this form.			
Last Name **		First Name **	
Address **		Telephone Number	
City **	State **	Zip Code **	Country of Origin / Bi-National ID #
DOB (mm/dd/yyyy) **	Sex **	Pregnant?	
		<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Unknown	
Race:	<input type="checkbox"/> White	<input type="checkbox"/> Black or African American	<input type="checkbox"/> Hispanic
	<input type="checkbox"/> American Indian / Native Alaskan	<input type="checkbox"/> Asian	<input type="checkbox"/> Non-Hispanic
	<input type="checkbox"/> Native Hawaiian / Pacific Islander	<input type="checkbox"/> Other	<input type="checkbox"/> Unknown
Date of Collection ** (REQUIRED)	Time of Collection	<input type="checkbox"/> AM <input type="checkbox"/> PM	Collected By
Medical Record #/Allen #/CCL	CDC ID	Previous DSHS Specimen ID	HMO / Managed Care / Insurance Company Name *
ICD Diagnosis Code ** (1)	ICD Diagnosis Code ** (2)	ICD Diagnosis Code ** (3)	Address *
Date of Onset	Diagnosis / Symptoms	Risk	City *
<input type="checkbox"/> Inpatient <input type="checkbox"/> Outpatient <input type="checkbox"/> Outbreak association		<input type="checkbox"/> Symptomatic	State *
Section 3. SPECIMEN SOURCE (TYPE) (** REQUIRED)		Zip Code *	
<input type="checkbox"/> Blood <input type="checkbox"/> CSF <input type="checkbox"/> Serum <input type="checkbox"/> Other: _____		Responsible Party (Last Name, First Name) *	
Recruitment storage conditions, date and time specimens were removed from storage: <input type="checkbox"/> FREEZER <input type="checkbox"/> DATE (mm/dd/yyyy) _____ <input type="checkbox"/> REFRIGERATOR <input type="checkbox"/> TIME (AM/PM) _____ <input type="checkbox"/> AM <input type="checkbox"/> PM		Insurance Phone Number * <input type="checkbox"/> Responsible Party's Insurance ID Number *	
Section 4. HIV/STO TESTING		*1 hereby authorize the release of information related to the services described here and hereby assign any benefits to which I am entitled to the Texas Department of State Health Services, Laboratory Services Section. Signature * _____ Date * _____	
<input type="checkbox"/> HIV Screen <input type="checkbox"/> Syphilis RPR Only (Confirmation Required) <input type="checkbox"/> Syphilis Screen ** (Confirmation by TP-PA: Justification Required)			
Section 5. HEPATITIS TESTING		Section 6. SEROLOGICAL REFERENCE TESTING	
<input type="checkbox"/> Hepatitis A IgM <input type="checkbox"/> Hepatitis A Total (IgM/IgG) <input type="checkbox"/> Hepatitis B Core Antibody IgM <input type="checkbox"/> Hepatitis B Core Total Antibodies (IgM/IgG) <input type="checkbox"/> Hepatitis B Surface Antibody <input type="checkbox"/> Hepatitis B Surface Antigen <input type="checkbox"/> Hepatitis C Antibody		<input type="checkbox"/> Brucella IgG <input type="checkbox"/> Q-Fever IgG <input type="checkbox"/> Ehrlichia IgG <input type="checkbox"/> Rocky Mountain Spotted Fever & Typhus Fever Panel IgG <input type="checkbox"/> Hantavirus IgM & IgG <input type="checkbox"/> Rubella IgM <input type="checkbox"/> Measles IgM <input type="checkbox"/> Rubella IgG <input type="checkbox"/> Measles IgG <input type="checkbox"/> Schistosoma IgG <input type="checkbox"/> Mumps IgG <input type="checkbox"/> Streptococcus IgG <input type="checkbox"/> Plague IgG <input type="checkbox"/> Tularemia IgG	
Section 9. CDC REFERENCE TESTS		Provide patient history on reverse side of form or attach to avoid delay of specimen processing <input type="checkbox"/> Chagas Disease <input type="checkbox"/> Leptospirosis <input type="checkbox"/> Cysticercosis <input type="checkbox"/> Paragonimiasis <input type="checkbox"/> Echinococcosis <input type="checkbox"/> VRDL (CSF only) <input type="checkbox"/> Fascioliasis <input type="checkbox"/> Other: _____ <input type="checkbox"/> HTLV-1	
FOR LABORATORY USE ONLY		Specimen Received: <input type="checkbox"/> Room Temp. <input type="checkbox"/> Cold <input type="checkbox"/> Frozen	

Laboratory Services Section: 1100 W 49th St Austin, Tx 78756

- Use the **G-2A Serology Specimen Submission Form** for Hepatitis B, C and HIV serology testing for patients in the TB program.
- **Section 1.** Ensure all information is updated and current.
- **Section 2.** Fill out completely. Patient name must match exactly the patient name labeled on the specimen container.
- **Section 3.** Specimen source must be provided.
- **Sections 4 and 5.** Requested test must be selected, or specimen is unsatisfactory for testing.
- **Section 7 and 8.** Fill out completely. Ensure *TB Elimination* is marked as "Payor Source" in Section 8.

Visit laboratory website for the most recent requisition version
https://www.dshs.texas.gov/lab/mrs_forms.shtm

F40-TB Elimination Specimen Submission Form - DSHS South Texas Laboratory (STL)

TEXAS South Texas Laboratory		Texas Department of State Health Services		CLIA #45D0503753 CAP #2148801		***DSHS LAB USE ONLY**	
P: (956) 364-8746 FAX: (956) 412-8794 https://www.dshs.texas.gov/lab/stl_lab.htm							
Section 1. SUBMITTER INFORMATION - (** REQUIRED)				Section 3. ORDERING PHYSICIAN INFORMATION - (** REQUIRED)			
Submitter/PI Number		Submitter Name		Ordering Physician's NPI Number		Ordering Physician's Name	
NPI Number		Address					
City		State		Zip Code			
Phone		Contact					
Fax		Clinic Code					
Section 2. PATIENT INFORMATION - (** REQUIRED)				Section 4. PAYOR SOURCE - (** REQUIRED)			
NOTE: Patient name on specimen is REQUIRED & MUST match name on this form & Medicare/Medicaid card.				1. Please do not use this form if not funded by the TB Elimination Program; use the F40-A specimen submission form.			
Last Name		First Name		MI		2. If the patient does not meet program eligibility (see form) for the test requested and no third party payor will cover the testing, the patient will be billed.	
Address		Telephone Number				3. Medicare generally does not pay for screening tests. Please refer to applicable Third party payer guidelines for instructions regarding coverage of tests, benefit limitations, medical necessity designations and/or prior authorization/Referral Notice (ARN) requirements.	
City		State		Zip Code		4. If Medicare or Medicare is indicated, the Medicare Medicare number is required. Please write it in the space provided below.	
DOB (mm/dd/yyyy)		Sex		SSN		5. If private insurance is indicated, the requesting information below is designated with an asterisk (*).	
Race		Ethnicity		Pregnant?		6. Check only one box to indicate whether we should bill the submitter: Medicaid, Medicare, Private Insurance or DSHS Program.	
<input type="checkbox"/> White <input type="checkbox"/> American Indian / Native Alaskan <input type="checkbox"/> Native Hawaiian / Pacific Islander		<input type="checkbox"/> Black or African American <input type="checkbox"/> Asian <input type="checkbox"/> Other		<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Unknown		<input type="checkbox"/> TB Elimination (1619)	
Date of Collection (** REQUIRED)		Time of Collection (**)		Collected By			
Medical Record Number		Area # / CD / CDC ID		Previous DSHS Specimen Lab Number		HMO / Managed Care / Insurance Company Name *	
ICD Diagnosis Code ** (I)		ICD Diagnosis Code ** (Z)		ICD Diagnosis Code ** (E)		Address *	
<input type="checkbox"/> Inpatient <input type="checkbox"/> Outpatient <input type="checkbox"/> Outbreak association		Diagnosis / Symptoms		Risk		City * State * Zip Code *	
Date of Onset (mm/dd/yyyy)						Responsible Party (Last Name, First Name) *	
Section 5. CHEM PANELS		Section 6. CHEMISTRY		Insurance Phone Number *		Responsible Party's Insurance ID Number *	
<input type="checkbox"/> Basic Metabolic Panel # (Sodium, Potassium, Chloride, CO2, Glucose, BUN, Creatinine, Calcium)		<input type="checkbox"/> Albumin <input type="checkbox"/> Alkaline Phosphatase <input type="checkbox"/> ALT (SGPT)		Group Name		Group Number	
<input type="checkbox"/> Comp Metabolic Panel # (Sodium, Potassium, Chloride, CO2, Glucose, BUN, Creatinine, ALT, AST, AB, Phos, TBL, AB, Total Protein, Calcium)		<input type="checkbox"/> AST (SGOT) <input type="checkbox"/> Bilirubin, Total <input type="checkbox"/> Blood Urea Nitrogen (BUN)		<input type="checkbox"/> Hemoglobin A1C <input type="checkbox"/> Magnesium <input type="checkbox"/> Protein, Total <input type="checkbox"/> Uric Acid		<input type="checkbox"/> Hepatic Function Panel (AB, ALT, AST, AB, Phos, TBL, DBI, Total Protein)	
<input type="checkbox"/> Renal Function Panel (Sodium, Potassium, Chloride, CO2, Glucose, BUN, Creatinine, AB, Calcium, Hemoglobin)		<input type="checkbox"/> Glucose <input type="checkbox"/> Hemoglobin A1C <input type="checkbox"/> Magnesium <input type="checkbox"/> Protein, Total <input type="checkbox"/> Uric Acid		<input type="checkbox"/> Thyroid stimulating hormone (TSH) <input type="checkbox"/> Thyroxine (T4), Total		<input type="checkbox"/> TB Panel: (ALT, AST, AB, Phos, TBL, BUN, Chol, Creatinine, Uric Acid)	
Section 8. HEMATOLOGY		Section 9. SPECIAL CHEMISTRY		<input type="checkbox"/> CBC automated with differential		<input type="checkbox"/> Thyroid stimulating hormone (TSH) <input type="checkbox"/> Thyroxine (T4), Total	
<input type="checkbox"/> CBC automated with differential				<input type="checkbox"/> Thyroid stimulating hormone (TSH) <input type="checkbox"/> Thyroxine (T4), Total		<input type="checkbox"/> TB Elimination (1619)	
NOTES: * = Fasting preferred for test. ** Document type & date specimens were removed from FREEZER/REFRIGERATOR in the lower right-hand box		FOR LABORATORY USE ONLY		Indicate removal from:		DATE TIME	
				<input type="checkbox"/> FREEZER <input type="checkbox"/> REFRIGERATOR		<input type="checkbox"/> Room Temp <input type="checkbox"/> Cold <input type="checkbox"/> Frozen	
				Specimen Received:			

Laboratory Services Section/South Texas Lab: 1301 S Rangerville Rd Harlingen, Tx 78552

- Use the **F40-TB Elimination Specimen Submission Form** for chemistry and hematology clinical blood samples.
- **Section 1.** Ensure all information is updated and current.
- **Section 2.** Patient name must match exactly the name labeled on the specimen container. Ensure date and time of collection are included.
- **Sections 3 and 4.** Fill out completely. Ensure *TB Elimination* is marked in Section 4 "Payor Source."
- **Sections 5, 6, 8 and 9.** Select each test requested. Magnesium may be ordered for patients on Bedaquiline only.

Visit laboratory website for the most recent requisition version
<https://dshs.texas.gov/lab/stlForms.htm>

F40-B Specimen Submission Form - DSHS South Texas Laboratory (STL)

F40-B Specimen Submission Form (Jan 2020)		Place DSHS Bar Code Label / Address-O-Graph Here													
Texas Department of State Health Services CLIA #45D0503753 CAP #2148801 www.dshs.texas.gov/lab/sto_tx_lab															
Section 1. SUBMITTER INFORMATION - (** REQUIRED) P: (956) 364-8746 FAX: (956) 412-8734															
Submitter/TPN Number: _____ Submitter Name: _____ NPI Number: _____ Address: _____ City: _____ State: _____ Zip Code: _____ Phone: _____ Contact: _____ Fax: _____ Clinic Code: _____		Section 3. ORDERING PHYSICIAN INFORMATION - (** REQUIRED) Ordering Physician's NPI Number: _____ Ordering Physician's Name: _____													
Section 2. PATIENT INFORMATION - (** REQUIRED) NOTE: Patient name on specimen MUST match name on this form & Medication/Medicaid card. Specimen must have two (2) identifiers that match this form. Last Name: _____ First Name: _____ MI: _____ Address: _____ Telephone Number: _____ City: _____ State: _____ Zip Code: _____ Country of Origin: _____ DOB (mm/dd/yyyy): _____ Sex: _____ SSN: _____ Pregnant? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Unknown Race: <input type="checkbox"/> White <input type="checkbox"/> Black or African American <input type="checkbox"/> Hispanic <input type="checkbox"/> Unknown <input type="checkbox"/> American Indian / Native Alaskan <input type="checkbox"/> Asian <input type="checkbox"/> Non-Hispanic <input type="checkbox"/> Unknown <input type="checkbox"/> Native Hawaiian / Pacific Islander <input type="checkbox"/> Other															
Section 4. PAYOR SOURCE - (** REQUIRED) 1. Reflex testing will be performed when necessary and the appropriate party will be billed. 2. If the patient does not meet program eligibility requirements for the test requested and no third party payor will cover the cost, the submitter will be billed. 3. Medicare generally does not pay for genetic testing. Refer to applicable Third party payor guidelines for instructions regarding coverage, benefits, limitations, medical necessity determinations, and Secondary Notice (ABN) requirements. 4. If Medicare or Medicare is indicated, the Medicare procedure number is required. Please write it in the space provided below. 5. If private insurance is indicated, the required billing information below is designated with an asterisk. 6. Check only one box below to indicate whether we should bill the submitter, Medicaid, Medicare, private insurance, or DSHS Program.															
Section 5. MEDICAL RECORD INFORMATION Medical Record #/Alem #/CUI: _____ CDC ID: _____ Previous DSHS Specimen Lab #: _____ ICD Diagnosis Code ** (1): _____ ICD Diagnosis Code ** (2): _____ ICD Diagnosis Code ** (3): _____ <input type="checkbox"/> Inpatient <input type="checkbox"/> Outpatient <input type="checkbox"/> Outbreak association: <input type="checkbox"/> _____		Medicaid (2) <input type="checkbox"/> Medicare (8) <input type="checkbox"/> Medicaid/Medicare #: _____ <input type="checkbox"/> IDEAS (3) <input type="checkbox"/> OPC <input type="checkbox"/> <input type="checkbox"/> RDS (1318) <input type="checkbox"/> Private Insurance (4) <input type="checkbox"/> <input type="checkbox"/> B/USST (1719) <input type="checkbox"/> TB Elimination (1619) <input type="checkbox"/> <input type="checkbox"/> HIV / STD (1608) <input type="checkbox"/> Other: _____													
Section 6. SPECIMEN SOURCE OR TYPE - (REQUIRED for Mycobacteriology specimens) <input type="checkbox"/> Abscess (site) <input type="checkbox"/> Lesion (site) <input type="checkbox"/> Throat swab <input type="checkbox"/> Blood <input type="checkbox"/> Lymph node (site) <input type="checkbox"/> Tissue (site) <input type="checkbox"/> Bone marrow <input type="checkbox"/> Nasopharyngeal <input type="checkbox"/> Urine <input type="checkbox"/> Bronchial washings <input type="checkbox"/> Plasma <input type="checkbox"/> Vaginal <input type="checkbox"/> CSF <input type="checkbox"/> Rectal swab <input type="checkbox"/> Wound (site) <input type="checkbox"/> Eye <input type="checkbox"/> Serum <input type="checkbox"/> Other: _____ <input type="checkbox"/> Feces/stool <input type="checkbox"/> Sputum <input type="checkbox"/> _____ <input type="checkbox"/> Gastric <input type="checkbox"/> _____															
Section 7. MYCOBACTERIOLOGY <input type="checkbox"/> AFB Culture <input type="checkbox"/> Identification of AFB isolate, DNA Probe <input type="checkbox"/> AFB Smear only <input type="checkbox"/> Identification, referred isolate, DNA Probe <input type="checkbox"/> AFB Concentration <input type="checkbox"/> MGIT Susceptibility (each drug) <input type="checkbox"/> Direct NAAT (M. tuberculosis) <input type="checkbox"/> MGIT Susceptibility (each drug) PZA <input type="checkbox"/> (Respiratory Diagnostic Specimen Only) <input type="checkbox"/> Conventional Susceptibility (each drug)															
Section 8. ZIKA, DENGUE, CHIKUNGUNYA <input type="checkbox"/> Zika, Dengue, and/or Chikungunya NCTE, Serology, PCR, or both will be performed at CHS and the testing methodology and specific viruses approved for testing will be based on clinical symptoms and epidemiological criteria. In some instances, specimens may also be forwarded to CDC for further testing.															
Section 9. LABORATORY TEST RESULTS SECTION - FOR LABORATORY USE ONLY <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>TEST</th> <th>NONREACTIVE</th> <th>REACTIVE</th> <th>TITER</th> </tr> </thead> <tbody> <tr> <td>RPR</td> <td></td> <td></td> <td></td> </tr> <tr> <td>TP-PA</td> <td></td> <td></td> <td></td> </tr> </tbody> </table> UNSATISFACTORY: <input type="checkbox"/> Broken in Mail <input type="checkbox"/> Lacked in Transit <input type="checkbox"/> No Specimen Received <input type="checkbox"/> Thawed <input type="checkbox"/> Hemolyzed <input type="checkbox"/> Name Discrepancy <input type="checkbox"/> Quantity Not Sufficient <input type="checkbox"/> Please resubmit: FOR LABORATORY USE ONLY: <input type="checkbox"/> Room Temp <input type="checkbox"/> Cold <input type="checkbox"/> Frozen				TEST	NONREACTIVE	REACTIVE	TITER	RPR				TP-PA			
TEST	NONREACTIVE	REACTIVE	TITER												
RPR															
TP-PA															

- Use the **F40-B Specimen Submission Form** for mycobacteriology and TB specimen testing.
- **Section 1.** Ensure all information is updated and current.
- **Section 2.** Patient name must match exactly the patient name labeled on the specimen container. Ensure date and time of collection are included.
- **Sections 3 and 4.** Fill out completely. Ensure *TB Elimination* is marked in Section 4 "Payor Source."
- **Section 5.** Specimen source must be provided. This will direct how the specimen is processed in the laboratory.
- **Section 7.** Requested test must be selected or specimen is unsatisfactory for testing.

Visit laboratory website for the most recent requisition version

<https://dshs.texas.gov/lab/stlForms.htm>