



Printing Newborn Screening Web Application (Neometrics) Labels with Firefox



Complete all patient demographic information fields

Once Patient Information has been completed select the 'Print Label' option.



When the Pop-Up window appears press CTRL+P to print.

From the Print Options select 'More Settings'

Uncheck the 'Print Headers and Footers' option

Under 'Margins' select 'Custom (inches)' and set all margins to 0.30

Print

Once Printed


The information will print on one label and Firefox will retain these settings as the default for future printing.

For the NBS Remote Data Website, visit:

<https://www.dshs.texas.gov/lab/remotedata.shtm>

For the Printer Configuration Manual, visit:

<https://www.dshs.state.tx.us/lab/nbsProviderWebInstallation.pdf>



1-888-963-7111 ext. 2484 | LabInfo@dshs.texas.gov | RemoteLabSupport@dshs.texas.gov