


1 Store vaccine at ideal temperature: 40°F




Never freeze refrigerated vaccine!

Exception: MMR can be stored in fridge or freezer



Report out of range temperatures immediately!

2 Record daily temperatures



Three Steps - Twice a Day: Temperatures should be checked and recorded first thing in the morning and before leaving at night.

- 1 Current Temperature:** The temperature that the refrigerator is right now.
- 2 Min/Max:** The coldest and warmest the refrigerator has been since you last reset the thermometer.
- 3 Reset:** The button you push after you have checked the Min/Max.



Best Practices

3 Take action if out of range!

- Contact your state or local health department immediately. Or if private vaccine call the manufacturer directly.
- Tell them the total amount of time the refrigerator was out of range.
- **Take your time** - Read and record temperatures accurately.
- **Make your mark!** Initial the log when recording temperatures.
- **Leave it blank** - if a temp was not recorded, leave the space blank!



U.S. Department of Health and Human Services
Centers for Disease Control and Prevention

Distributed by



Texas Department of State Health Services
Immunization Branch

Stock No. E11-14484
Rev. 10/2014


CS243541-A Revision Jan 24, 2014

Visit www.cdc.gov/vaccines/SandH
for more information, or your state health department.

Test Your Knowledge


Review the below temperature readings and select the correct answer.

1




- A. Current temp and min/max are within range - no action necessary
- B. Current temp is within range, min/max out of range - take action
- C. Current temp is within range, min/max out of range - no action necessary
- D. Current temp and min/max are out of range - take action

2




- A. Current temp and min/max are within range - no action necessary
- B. Current temp is within range, min/max out of range - take action
- C. Current temp is within range, min/max out of range - no action necessary
- D. Current temp and min/max are out of range - take action

3



- A. Current temp and min/max are within range - no action necessary
- B. Current temp is within range, min/max out of range - take action
- C. Current temp is within range, min/max out of range - no action necessary
- D. Current temp and min/max are out of range - take action

4



- A. Current temp and min/max are within range - no action necessary
- B. Current temp is within range, min/max out of range - take action
- C. Current temp is within range, min/max out of range - no action necessary
- D. Current temp and min/max are out of range - take action

5

Take action means (circle any that apply)

- A. Remove all vaccines that were out of range and discard them
- B. Call the state/local VFC program (or manufacturer if private vaccine) for guidance
- C. Notify practice vaccine coordinator to change the temperature controls to get the refrigerator back in range
- D. Thaw any vaccines that got frozen for 45 minutes